

Baird microwave ablation (MWA) generator

Integrated control, cooling, and monitoring for image-guided soft-tissue ablation.



Intuitive control with real-time feedback

Set power, time, and output mode while viewing pump status, needle-shaft temperature, and alerts on one screen.

No grounding pad required

MWA does not require a patient-return electrode, eliminating the grounding-pad placement step in procedure setup.¹

Integrated cooling with output interlock

Built-in pump control supports water or saline circulation through Baird MWA needles. Microwave output remains disabled until cooling circulation is running.

Baird MWA generator at a glance



Product specifications

Microwave frequency	2450 MHz ± 25 MHz
Adjustable output	0-70 W*
Output modes	Foot · Continuous · Pulse
Cooling	Built-in peristaltic pump; distilled water or saline
Needle-shaft temperature safeguard	Audible alert above 41 °C; automatic microwave output cutoff at 45 °C
Power requirements	100-240 V~, 50/60 Hz, 500 VA
Dimensions (L x W x H)	41 × 43 × 15.5 cm / 16.1 × 16.9 × 6.1 in
Weight	11 kg / 24.3 lb

*Maximum output varies by compatible Baird MWA needle series; T-series up to 60 W, L-series up to 70 W.

Manufactured by Betters (Suzhou) Medical Co., Ltd. U.S. sales and support by Baird Medical LLC.

Indications & safety

For healthcare professionals. Rx only. The Microwave Ablation System (BD-GT) is indicated for coagulation (ablation) of soft tissue and is not intended for cardiac use. Contraindicated in patients with massive bleeding tendency; not recommended for pregnant women; use caution in patients with severe coagulation dysfunction. Potential complications of thyroid ablation reported in published literature include pain, bleeding or hematoma, swelling, thyroid-function changes, voice changes, and injury to nearby structures.² See the current IFU for complete indications, contraindications, warnings, precautions, and operating instructions. Individual outcomes may vary.

Reference: 2. Sinclair CF, et al. Thyroid. 2023;33(10):1150-1170.



12 E 49th St., 18th floor
New York, NY 10017

www.bairdmed.com | info@bairdmed.com |
+1 (888) 508-6228



©2026 Baird Medical. All rights reserved.
For distribution in the USA only.